

Evaluation Report [*Final Amended*]

Name of the Test Evaluated	:	Colour Scale for Hemoglobin (Hemocheck)
Name of Manufacturer	:	Allied Health Sciences Pvt.Ltd, 101-C, Kundan Plaza, Ashram, New Delhi- 110014
Batch Number	:	18/2727
Expiry Date	:	December 2018
Evaluated During	:	February 2014
Evaluation Protocol	:	<ul style="list-style-type: none"> ➤ Reference Range of the test using 20 healthy individuals ➤ Precision <ul style="list-style-type: none"> ▪ Intra-Assay Variation (One specimen assayed 20 times in same run) ▪ Inter-Assay Variation (One Specimen assayed 5 times in 5 runs) ➤ Method Correlation with automated cell counter (LH-750) using 25 specimens

Note: This report is amended as there were typographical errors in the first report.

Results	:	Obtained Reference Range - All the 20 normals showed hemoglobin values within the recommended range (> 12.0 g/dL)
	:	Precision – Intra-Assay Variation – 0.00% and Inter-Assay Variation – 0.00%.
	:	Method Correlation - Regression analysis (r = 0.95), indicates acceptable correlation of results with automated Cell Counter LH750
Conclusion	:	Performance of Colour Scale test kit for Hemoglobin was found to be acceptable.

Prepared by

07/03/2014
Mr Mathew Muthukutty
Junior Research Officer

Rr 2/3/2014
Ms Rashmi Khadapkar
DGM R&D Operations

Reviewed by:

Maria Borges 07/03/2014
Dr Maria Borges
Chief Quality Officer

Approved by:

Dr. B.R. Das 07/13/14
Dr. B.R. Das
President - Research & Innovations

Report Number : R&D-HEM020v1

Page 2 of 2